Operating manual

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ACROBAT  
*Swing*

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**1 Instructions for safe usage ..... 4**

1.1 Explanation of symbols ..... 6

1.1.1 Symbols in the operating instructions ..... 6

1.1.2 Marking on the equipment ..... 6

1.1.3 Information on the rating plate (Fig.1) ..... 8

1.2 General Safety indications ..... 9

1.2.1 Standards and guidelines ..... 9

1.2.2 Intended purpose ..... 9

1.2.3 Incorrect use ..... 9

1.2.4 Contraindications ..... 9

1.2.5 Ambient conditions for operation and storage ..... 10

1.2.6 Cleaning and Disinfection ..... 10

**2 Using the ACROBAT Swing Equipment ..... 11**

2.1 Components of the equipment (Fig. 2) ..... 11

2.2 Usage (Fig. 2) ..... 11

2.3 Disposal ..... 14

**3 Cleaning and Disinfection ..... 15**

3.1 General Safety indications ..... 15

3.1 Cleaning ..... 15

3.2 Disinfection ..... 15

**4 Tests and Maintenance ..... 16**

4.1 Tests ..... 16

4.2 Maintenance ..... 16

**5 Maintenance ..... 17**

5.1 Adjusting the spring force (Fig. 3) ..... 17

5.2 Greasing the securing segment ..... 19

**6 Dismantling and mounting for service ..... 20**

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# 1 Instructions for safe usage

## Dear operator,

Please note that certain operations must only be carried out by personnel with the required qualifications:

- **The device may only be operated by persons with medical knowledge,**
- **The device may only be cleaned by trained cleaning personnel.**
- These operating instructions describe the ACROBAT Swing tripod-mounted, wall-mounted and ceiling-mounted device.  
The end devices (e.g. luminaires for medical diagnosis, lamps etc) are supplied by other manufacturers. For using such equipment, please refer to the operating instructions of those manufacturers.
- Please read these operating instructions carefully before using the equipment. You will then be able to exploit all the advantages that the equipment offers, and at the same time, protect yourself and others from injury.
- The device is only intended for use for the purposes described in the operating instructions. Any other application can cause injury or death and jeopardize the equipment and other material assets of the user.
- If problems occur that are not adequately described in these operating instructions, please do contact your supplier for your own safety.

### Hand-over declaration for the user:

I have read and understood the operating instructions. I am conversant with the safety instructions. I have been adequately instructed in preventive measures for avoiding the dangers and am following these operating rules when using the machine.

**Please keep these operating instructions carefully in the vicinity of the equipment so that you can later refer to safety instructions and important information at any time.**

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# 1 Instructions for safe usage

## Dear user:

Please note that certain work may only be carried out by suitably qualified personnel.

- **The work in chapters 5 and 6 may only be carried out by a hospital technician or a person with comparable qualifications.**
- This equipment is state-of-the-art and is operationally safe. Nonetheless, it may present some dangers, especially if it is operated by inadequately trained personnel or if it is used improperly or not used for its intended purpose.
- With the help of these operating instructions, please instruct the persons responsible for operating and cleaning in the working, operation and maintenance and care of this machine.
- The contents of these operating instructions can be changed at any time without prior notice.
- In case of translations in other languages, ONLY the German version of these instructions is legal and binding.
- Unauthorized changes made by the user or modifications to the equipment are not permitted for reasons of safety.

**Please keep these operating instructions carefully in the vicinity of the equipment so that you can later refer to safety instructions and important information at any time.**

The Acrobat Swing appliances must be installed prior to use. Installation instructions are included in the scope of delivery of the product.

## Trademarks:

- ACROBAT Swing® is a registered trademark of Ondal.
- All other trademarks mentioned in these operating instructions are the exclusive property of their respective owners and are acknowledged as such.

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# 1 Instructions for safe usage

## 1.1 Explanation of symbols

In the operating instructions, important instructions have been marked with symbols. The pictograms have the following meanings:

### 1.1.1 Symbols used in the operating instructions:

**DANGER!**  
Non-observance of these will result in serious or even fatal injuries.

**WARNING!**  
Non-observance of these may cause serious or fatal injuries.

**CAUTION!**  
Non-observance of these may cause medium to light injuries or material damage.

**IMPORTANT!**  
Provides usage tips and useful information.

 CE-Conformance mark

### 1.1.2 Graphic symbols on the appliance and/or on the packaging:

**Observe the Operating Instructions:**

Read these Operating Instructions carefully prior to initial operation of the pendant system. This ensures that you benefit from all the advantages of the pendant system and prevents any risk of injury or damage.

**Observe the maximum loading capacity:**

Warns of exceeding the maximum loading capacity (payload) on the pendant system of the adaption or the end device (e.g. flat screen, OR lamp, etc.).

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## 1.1 Explanation of symbols

 CE mark:

Ondal declares that the products comply with the relevant regulations set forth in the applicable European Directives.

 Recognised cUL component:

This component has been recognised by Underwriters Laboratories Inc. Representative samples of this product have been reviewed by UL and comply with the applicable requirements.

**Atmospheric pressure:**

indicates the permissible atmospheric pressure values in a range from 500 hPa to 1060 hPa for transport and storage.

**Relative humidity:**

indicates the permissible humidity values in a range from 10% to 75% for transport and storage.

**Ambient temperature:**

indicates the permissible ambient temperature values in a range from -25 °C to 70 °C for transport and storage.

**Risk of tilting:**

The appliance has a maximum loading capacity of 39 Nm (48 Nm special variant). If the maximum loading capacity is exceeded, there is a risk of the appliance tilting, hitting persons and causing severe injury.

- The maximum loading capacity of 39 Nm (48 Nm special-variant) must not be exceeded.
- Do not climb on the tripod and do not lean against it.
- Do not attach any additional loads to the top and side panels of the appliance.

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# 1 Instructions for safe usage

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## 1.1 Explanation of symbols

### Moving the appliance to a different location

Be aware that the spring arm must be placed to its lowest position (arrow 1) before moving the tripod version to a different location.

- Makes sure that the lockable rollers are unlocked.
- Do not pass over uneven flooring, thresholds, landings on lift doors and other obstacles.
- Make sure that you move the appliance at a reasonable speed in order to be able to stop it or move it around obstacles at any time.
- Be aware of sloping ground.

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# 1 Instructions for safe usage

## 1.1.3 Information on the rating plate Fig. 1

### The rating plate is attached to the following parts:

At the bottom of the tripod base. Under the tripod handle or at the top near the connection to the spring arm. At the top of the spring arm near the joint. At the top of the wall panel. On the extension arm near the connection to the wall bearing or ceiling mount fixture. On the flange plate under the canopy and/or on the ceiling tube near the connection to the extension arm.

### Serial number

- The rating plate indicates the serial number (SN) and the reference number (REF).

### Power supply

- The rating plate indicates voltage and current data (see also Technical Data).


### Load bearing capacity / weight

- For example, 30 Nm indicates the maximum permissible carrying load on the tripod. 12 kg indicates the weight of the tripod base (see also Technical Data).

### Date of manufacture

- The digits 1 to 4 of the serial number (SN) indicate the date of manufacture of the pendant system.
- The first two digits indicate the week of manufacture, e.g. 14 = calendar week 14.
- The following two digits indicate the year of manufacture, e.g. 11 = 2011.
- The letter in the 5th position indicates the factory, e.g. H = Hünfeld.
- The digits following the letter indicate the serial number. The information and illustrations serve as examples.

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
# 1 Instructions for safe usage

## 1.2 General safety instructions

The specific safety instructions in the following chapters must be adhered to.

### 1.2.1 Standards and guidelines

The appliance complies with the safety requirements of the following standards and directives:

- Medical Devices Act (Medizinproduktegesetz - MPG);
- MDD 93/42/EEC, 2007 – Medical Device Directive;
- EN 60601-1: 2006 - Medical Electrical Equipment - Part 1: General Requirements for Basic Safety and Essential Performance.

### 1.2.2 Intended purpose:

- The Acrobat appliances have been designed for the carrying and positioning of medical examination lamps and flat screens and also for the power supply of the devices. Depending on the transformer version, e.g. 13,2V, 13,5V, 14V, 24V (without transformer, 230V), 50/60Hz are available (not optionally adjustable).
- The appliance is suitable for continuous operation.
- The maximum load indicated on the rating plate of the individual appliance version must not be exceeded (see Chapter 6.2 "Mounting the end device").

### 1.2.3 Incorrect use

- The maximum loading capacity of the appliance and its components as specified in "Chapter 7, Technical Data" must not be exceeded.

### 1.2.4 Contraindications

- The ACROBAT Swing pendant system must not be used close to strong magnetic fields.
- No BF or CF application parts in accordance with IEC 60601-1 may be directly connected to the pendant system.

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# 1 Instructions for safe usage

## 1.2.5 Ambient conditions for operation and storage

Different conditions are required for the operation and interim storage of the appliance

### Ambient conditions for storage and transport

The following storage conditions apply for storage times of up to 15 weeks:

- Ambient temperature: -25 °C to 70 °C;
- Relative humidity: 10 % to 75 %;
- Atmospheric pressure: 500 hPa to 1060 hPa.

Store only in indoor rooms; after this time, the values specified for the ambient conditions for operation apply.

### Ambient conditions for operation

- Ambient temperature: 10 °C to 40 °C;
- Relative humidity: 30 % to 75 %;
- Atmospheric pressure: 700 hPa to 1060 hPa.

## 1.2.6 Cleaning and disinfection

### Cleaning

 **WARNING - Risk of contamination and infection of the patient**

Parts of the pendant system and the adaptations are made of plastic; solvents can dissolve plastic materials. Strong acids, bases and agents with an alcoholic strength of more than 60% can lead to the plastic materials becoming brittle.

Damaged parts can fall into open wounds. Liquid cleaning agents penetrating the pendant system and the adaptations can drip into open wounds.

### Disinfection

 **WARNING - Health hazard**

Disinfectants can contain substances hazardous to health which, when in contact with the skin and eyes, can cause injuries or affect the respiratory organs when inhaled. Observe the protective measures:

- Observe the hygiene regulations.
- Adhere to the disinfectant manufacturer's instructions.
- Perform surface disinfection every working day and in case of contamination.

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## 2 Using the ACROBAT Swing devices

 **NOTE: Qualifications of the personal:**  
The device must be operated by medically qualified expert technicians in hospitals and in physician's clinics.

### 2.1 Components of the devices Fig. 2

Ceiling device:	Wall-mounted device:	Tripod device
<ul style="list-style-type: none"><li>• Ceiling bearing (1)</li><li>• Space tube (2)</li><li>• Extension arm (3)</li><li>• Spring arm (4)</li></ul>	<ul style="list-style-type: none"><li>• Power Plug (5)</li><li>• Wall bearing (6)</li><li>• Extension arm (7)</li><li>• Spring arm (8)</li></ul>	<ul style="list-style-type: none"><li>• Spring arm (9)</li><li>• Tripod tube (10)</li><li>• Tripod base (11)</li><li>• Rollers (12)</li><li>• Power plug (13)</li></ul>

### 2.2 Usage Fig. 2

 **NOTE:**  
Because of the danger of the tripod mounting tipping over, care must be taken that it should not roll over any objects on the floor or the connecting cable.

1. Move the tripod device to the location of usage and **secure the two fastenable rollers (12)**.

#### Mounting / dismantling:

 **CAUTION - Electrical shock:**  
**To prevent the risk of electric shock, the appliance must be connected to a supply network with a protective conductor.**

- The appliance must be connected in such a way that it can be disconnected from the mains at all poles and at the same time.

 **WARNING - Electrical shock:**  
**In case of contact with damaged electrical parts, there is a danger of an electrical shock. In case of any damage to the power supply cable or plug, do not connect the equipment to the power supply and immediately inform the service department.**

2. Check the **power plug and the connection cable (5 / 13)** for damage.
3. To prevent the risk of electric shock, the mains plug may only be plugged into a socket outlet with grounding contact (Schuko socket) which has been correctly installed and earthed.

**Combination with other medical products:**  
The spring arm is combined with end devices of third-party

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## 2 Using the ACROBAT Swing Equipmente


manufacturers. Read the Operating Instructions provided by the third-party manufacturer and in particular the relevant pages with information on the operation of the end device.

The spring arm is only intended for use with Ondal extension arms, wall bearings and ceiling mount fixtures. If third-party products are combined with Ondal products, the manufacturer/marketer of the ME system must provide a declaration in accordance with Article 12 of 93/42/EEC (Medical Device Directive, MDD).

For technical assistance for creating a medical-electrical (ME) system refer to DIN EN 60601-1:2006, Part 16.

**Power packs intended for the supply of end devices must ensure electrical isolation and provide 2 protective measures in accordance with EN 60601-1.**

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
## 2 Using the ACROBAT Swing Equipmente

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### Patient Environment

Medical electrical devices which contain tangible parts and are attached to the tripod and positioned within the patient environment must provide two Means of Patient Protection (MOPP) in accordance with EN 60601-1.

If tangible parts are positioned outside the patient environment, two Means Of Operator Protection (MOOP) in accordance with EN 60601-1 must be provided.

The dimensions in the figure show minimum extent of the patient environment in a free surrounding.

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
Figure / Abildung 3

## 2 Using the ACROBAT Swing Equipment

 **NOTE: Qualifications of the personal:**  
The equipment must be operated by medically qualified expert technicians in hospitals and in physician's clinics.

### 2.3 Disposal

 **WARNING - Sudden release of spring tension**  
**The spring arm is equipped with a compression spring. When dismantling the spring arm, the compression spring suddenly releases its tension and can lead to severe injury:**

- Do not dismantle the spring arm for disposal.

#### RoHS conformity

- The appliance complies with the requirements of the EU Directive 2002/95/EC ("RoHS Directive" on the restricted use of certain hazardous substances in electrical and electronic appliances).
- To prevent environmental damage and personal injury, we therefore request you to contact us or your authorised service partner if you intend to take the appliance out of operation for the purpose of disposal.
- The appliance must be disposed of at a suitable collection point for the recycling of electrical and electronic devices in accordance with country-specific regulations.

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### 3.1 General safety instructions

#### **WARNING - Electrical shock:**

The appliances can carry an electric current and must be treated with the utmost care during cleaning and disinfection.

- If a mains plug exists, pull the mains plug.
- Do not apply spray cleaning and/or spray disinfection.
- Do not spray liquid into power sockets, gas sockets or appliance openings and prevent the penetration of liquids.

### 3.2 Cleaning

#### **Follow the safety instructions**

1. Follow the general safety instructions prescribed in "Chapter 4.1."

#### **Recommended cleaning agents**

Use a mild soap solution or a regular dishwashing product..

2. Wipe the surfaces of the appliances with a moderately moist cloth; add a mild soap solution (dishwashing product) if required.
3. Afterwards, carefully wipe the surfaces dry with a clean cloth.

### 3.3 Disinfection

#### **Follow the safety instructions**

Follow the general safety instructions prescribed in "Chapter 4.1."

#### **Disinfection method**

Wiping disinfection is the standardised disinfection method prescribed for the pendant system. Hygiene regulations and related safety instructions for the disinfection methods to be applied must be defined by the operator.

- In case of contamination with potentially infectious material (e.g. blood, body secretion or excrement) the surfaces must be immediately and specifically disinfected.
- Make sure that you apply the disinfectant in the correct concentration.
- For surface disinfection do not spray, but wipe, the surfaces.
- Wiped surfaced may only be used after the disinfectant has dried.

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## 4.1 Tests

### NOTES - Tests:

The user must check all ACROBAT Swing products with regard to the following points:

### WARNING - Electrical shock:

For all testing work, power off the equipment, pull out the plug and secure the equipment from being switched on again accidentally.

### Repeated inspections:

Repeated inspections must be carried out in accordance with DIN EN 62353.

### Twice a year:

- Damage to paint
- Cracks in plastic parts
- Deformation of the support system

### Annually:

- Check the fastening bolt on the lower side of the tripod base and if required, tighten it.
- Check the securing segment according to chapter 5.2, „Greasing the securing segment“, and grease it.

In case of any damage or faults, please inform your supplier.

## 4.2 Maintenance

### NOTE - annual maintenance:

All ACROBAT Swing products should be checked once a year with regard to the following, and maintenance carried out accordingly:

- Function test
- Electrical safety test
- Checking the support system

Your supplier is informed and trained as regards the scope of the maintenance work.

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# 5 Maintenance

 **NOTE - Qualification of the personnel:**  
The maintenance and care work must be carried out by a hospital technician (or a comparably qualified person) only.

## 5.1 Adjusting the spring force Fig. 4

Like every technical component, springs are subjected to natural wear. Thus, the spring force may give way and reduce after long operation, and has to be re-adjusted.

**Adjust the spring force in such a way that, the spring arm with end device comes to a standstill in every desired position.**

 **CAUTION - Destruction of the spring arm:**  
The spring tension is adjusted in the upper end position.

1. Remove the left-hand **joint cover half (2)** from the spring arm (when looking from the end device). To do this, gently push the joint cover half out of the groove in the spring arm joint using a small slotted screwdriver.
2. Move the end device to the upper end position.
3. Insert the slotted screwdriver into the **drill hole (1)** and adjust the spring tension.
4. Mount the joint cover and make sure that it snaps into place.

 **CAUTION - Destruction of the spring arm:**  
Driving in the brake screw too deeply destroys the spring arm. Tighten the brake screw with care while continually checking the braking force.

**If the spring arm drops - the spring force is too low.**

- The adjusting screw must be rotated to the left (in the counter-clockwise direction).

**If the spring arm goes upward - the spring force is too high:**

- the adjusting screw must be rotated to the right (in the clockwise direction).

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
Figure / Abildung 4

## 5.2 Greasing the securing segment

1. Dismantle the end-device as described in chapter 6.1.
2. Check the securing segment for a minimum thickness of 1.5 mm and if required, replace it.
3. Grease the securing segment and the pin of the end-device.
4. Mount the end-device as described in chapter 6.2, „Mount end-device“.

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## 6 Dismantling and mounting for servicing

**NOTE: Qualifications of the personnel**  
The assembly / dismantling must be carried out by a hospital technician (or a comparably qualified person).


### 6.1 Dismantling the end-device Fig. 5

#### **WARNING - danger of injury**

The spring arm is under a high spring tension. If the end device is not dismantled in the uppermost position of the spring arm, it rocks upwards and can result in serious injury. Dismount the end device only when the spring arm is in the uppermost position.

1. Power off the equipment/ pull out the power plug and secure it from being switched on.
2. Bring the spring arm in the uppermost position.
3. Rotate the **brake screw (1)** so that the **sleeve (2)** can be rotated.
4. Rotate the **sleeve (2)** through 180 degrees, till the **securing segment (3)** appears in the **slot (4)**.
5. Hold the **end device (5)** firmly and take out the **securing segment (3)** with a narrow slotted-screw screwdriver. Disconnect the **end device (5)**.
6. Put on the protective cap.

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
## 6 Dismantling and mounting for servicing

 **NOTE - Qualifications of the personnel:**  
The dismantling / assembly must be carried out by a hospital technician (or a comparably qualified person).

### 6.2 Mounting the end-device Fig. 6

 **NOTE:**  
For the tripod version, the load torque on the spring arm must not exceed 39 Nm (48 Nm special variant); for the wall-mounted and ceiling-mounted versions, it must not exceed 30 Nm (57 Nm special variant).

1. Power off the equipment / pull out the power plug and secure it from being switched on again.

 **WARNING - Danger of injury:**  
The spring arm, when it is pressed downward, can jump up suddenly and cause injury.  
Nobody should be present in the swiveling range of the spring arm during the installation of the end-device.

 **CAUTION - Damage to the device:**  
The spring arm will be damaged if it is pulled beyond the lower / upper stop. Grab the spring arm at the front end (not at the rear end) and slowly pull it downward until it reaches the stop.

2. Grab the spring arm at the front end and slowly pull it downward until it reaches the stop. Take out the protective cap from the spring arm opening.
3. Push on the **sleeve (1)** onto the spring arm in such a way that the **slots (2)** are covered.
4. Push in the **end-device (3)** and the **securing segment (4)** completely in the **slot (2)**, so that the securing segment is guided in the **groove (5)**.
5. Rotate the **plastic sleeve (1)** through 180 degrees, and tighten the **braking screw (6)** till the **end-device (3)** is braked.
6. Check the secure seating of the **end-device (3)**.

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
# 6 Dismantling and mounting for servicing

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Figure / Abildung 6

## 6 Dismantling and mounting for servicing

Acrobat Swing appliances are only intended for the connection of medical examination lamps, OR lamps and flat screens bearing the CE mark and - depending on the individual version - for the load torque indicated on the rating plate.

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The manufacturer / marketer of the ME system must provide a declaration in accordance with Article 12 of 93/42/EEC (Medical Device Directive, MDD).

 **NOTE - Built-in overload protection for versions with toroidal transformer:**

To protect the transformer and the end device, the fuses are installed.

### 6.3 Change the fuses Fig. 7

 **Warning - Electrical shock:**

For all maintenance work, switch off the equipment, pull out the main plug and secure it from being switched on again.

The technical specification for the fuses can be found on page 26 of this document and on the label inside the cover.

 **Caution - Damage to the Equipment**

Apply only the specified fuses.

Change the fuses corresponding to the following steps:

**Open the cover**

1. Loosening the **screw (3)**, don't unscrew completely.
2. Push up **clamp ring (2)** and **cover (1)** and secure them.

**Change fuse**

3. Change broken **fuse (4)**.

**Close the cover**

4. Push down **clamp ring (2)** and **cover (1)**.
5. Tighten **screw (3)**.

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



Figure / Abildung 7

# 7 Technical Data

Electrical data	Ceiling-unit	Wall-unit	Stativ-unit
Rated voltage, depending on model	120 V / 230 V	120 V / 230 V	120 V / 230 V
Maximum power*	90 W	80 W	90 W
Nominal frequency	50 / 60 Hz	50 / 60 Hz	50 / 60 Hz
Output current at 230 V *	bis zu 10 A	bis zu 10 A	bis zu 10 A
Typ of fuses: 230 V: T 10 L 250 V prim. 115V/sec. 14V/60VA: T1,0A/250V prim. 230V/sec. 14V/60VA: T0,5A/250V prim. 103V-127V/sec. 13,2V-14V/60VA: T1,25A/250V prim. 213V-247V/sec. 13,2V-14V/60VA: T0,63A/250V prim. 230V/sec. 13,5V/60VA: T0,5A/250V prim. 115V/sec. 24V/90VA: T1,6A/250V prim. 230V/sec. 24V/90VA: T0,8A/250V			
* depending on version			
Weight in kg	Ceiling-unit	Wall-unit	Stativ-unit
Spring arm (standard)	1,8 kg	1,8 kg	1,8 kg
Extension arm	1,9 kg	1,9 kg	—
Ceiling mounting	3,0 kg - 7,9 kg	—	—
Wall bearing	—	1,8 kg	—
Stativ			3,5 kg
Stativ base			6,5 kg
Stativ base (weighted)			12,0 kg 15,0 kg 18,0 kg
Classification according to EN 60601-1			
Device of protection class I			
IP-Classification according to IEC-60529			IP20
Suitable for continuous operation			
Classification according to guidelines 93/42/EWG, appendix IX			
Class I			
EU-Conformance			
The Acrobat Swing equipment conforms to the specifications of the guideline 93/42/EWG (medical products guideline). The Acrobat Swing equipment is only intended for connecting medical examination lamps, operation lamps and flat-screens with CE-markings and - je nach Ausführung - a maximum of 57 Nm self-weight. The person operating the equipment must submit a declaration as per article 12 of guideline 93/42/EWG (medical products guidelines). Another conformance evaluation must be carried out for other end devices.			

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### Important Information

- Inspections must be performed by trained service personnel.
- The inspection intervals must be observed.
- This inspection plan is only valid when combined with the Ondal Installation and Operating Instructions which must be complied with complementary to the inspections.
- After 10 years, the functional inspection must be performed once a year.

**At the specified intervals, the pendant system must be inspected for the following points by Ondal, a company authorised by Ondal, or by personnel with the corresponding qualification:**

### Visual inspection (once a year)

- The parts of the pendant system are not deformed\*
- The system is free from defects in paint wor
- The plastic parts are available and in position\*
- The plastic parts are free from cracks\*
- All rating plates are in position and legible

Period of use in years									
1	2	3	4	5	6	7	8	9	10
OK	NOK	OK	NOK	OK	NOK	OK	NOK	OK	NOK
<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

### Functional inspection (every 2 years)

- Free rotatability/limit stops (depending on the version)\*\*
- horiz. a. vert. joints smooth running, grease if required
- Vertical lift correct, readjust if required\*\*
- Check and grease securing segment\*
- Securing ring in position and shape (extension/spring arm)
- Load comp./spring tension correct, readjust if required
- Collision damage
  - all welding points free from cracks\*\*
- Inspection of protective conductor transition resistance\*\* (only applicable if current conducting cables are integrated)
- After 10 years, the cables and hoses (if existing) must be inspected and replaced if required\*\*

2	4	6	8	10	
OK	NOK	OK	NOK	OK	NOK
<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

### Confirmation of the inspection performed

The inspection tasks mentioned above were executed incl. the necessary adjustments and safety inspection:

1st year		6th year	
_____	Signature/stamp	_____	Signature/stamp
Date		Date	
2nd year		7th year	
_____	Signature/stamp	_____	Signature/stamp
Date		Date	
3rd year		8th year	
_____	Signature/stamp	_____	Signature/stamp
Date		Date	
4th year		9th year	
_____	Signature/stamp	_____	Signature/stamp
Date		Date	
5th year		10th year	
_____	Signature/stamp	_____	Signature/stamp
Date		Date	
_____	Signature/stamp	_____	Signature/stamp

\* Damaged or deformed components should be replaced as a precaution. Please contact the spring arm supplier.  
 \*\* If one of these problems is detected during the inspection, the spring arm must be taken out of service immediately as a precaution in order to prevent any further damage to persons or equipment. Inform the system suppliers immediately.  
 The medical device handbook which belongs to each medical device and is prescribed in accordance with the German Medical Product User Regulations (MPBetreibV) must be available on site. Service and maintenance work and also safety inspections shall be documented in this medical device handbook. Inspection reports, including the present report, must be filed in the corresponding medical device handbook.

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- Acrobat Swing     Acrobat 2000     Acrobat 3000  
 Acrobat 77        Acrobat LCH        OndaSpace

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© Ondal Medical Systems GmbH, 2012