



Hillrom™

Welch Allyn® Multi-Parameter Vital Signs Devices



Spot Vital Signs® 4400 Device



Connex® Spot® Monitor



Connex Vital Signs Monitor



Connex Integrated Wall System

WORKFLOWS	Spot	Spot, Office, Intervals	Spot, Office, Intervals, Continuous	Spot, Office, Intervals, Continuous
DEVICE TYPE	Basic	Intermediate	Advanced	Advanced
PATIENT POPULATIONS	Adult, Pediatric	Adult, Pediatric, Neonatal*	Adult, Pediatric, Neonatal*	Adult, Pediatric, Neonatal*
DEVICE ALARMS	—	Yes	Yes	Yes
THIRD-PARTY ALARM NOTIFICATION	—	—	Yes	Yes
PARAMETERS	Blood Pressure, Pulse Rate, Temperature, SpO ₂	Blood Pressure, Pulse Rate, Temperature, SpO ₂ , RRp	Blood Pressure, Pulse Rate, Temperature, SpO ₂ , etCO ₂ , SpHb, ECG, RRa	Blood Pressure, Temperature, SpO ₂ , SpHb, ECG, Otoscopy, Ophthalmoscopy
SUREBP® TECHNOLOGY (15-SEC)	Yes	Yes	Yes	Yes
BLOOD PRESSURE AVERAGING	Yes	Yes	Yes	Yes
DIAGNOSTIC HANDLES	—	—	—	Otoscope, Ophthalmoscope
SPO ₂ OPTIONS	Nonin®	Masimo®, Nellcor®, Nonin	Masimo, Nellcor	Masimo, Nellcor
TEMPERATURE OPTIONS	SureTemp® Plus	SureTemp Plus, Braun PRO 6000	SureTemp Plus, Braun PRO 6000	SureTemp Plus, Braun PRO 6000
EARLY WARNING SCORING	—	Up to 6 protocols	Up to 6 protocols	Up to 6 protocols
CUSTOM DOCUMENTATION FIELDS	5 Fixed and Selectable	Up to 20	Up to 20	Up to 20
WIRELESS CONNECTIVITY	—	WiFi®, Bluetooth®	WiFi	WiFi
WIRED CONNECTIVITY	USB	USB, Ethernet	USB, Ethernet	USB, Ethernet
DISPLAY	7" Color LCD, Touchscreen	7" Color LCD, Touchscreen	8.9" Color LCD, Touchscreen	8.9" Color LCD, Touchscreen
SCALE INTEGRATION	—	Yes	Yes	Yes
MOUNTING OPTIONS	Desktop Stand, Classic Stand	Power Management Stand, Classic Stand, Work Surface Stand, Wall Mount, Integrated Wall System, Desktop Stand	Cable Management Stand, Classic Stand, Wall Mount	Wall Mounted Unit
ACCESSORY OPTIONS	—	Barcode, RFID	Barcode, RFID	Barcode, RFID
MEMORY (# READINGS)	50	400	400	400

*Not all parameters may be used with neonates, plus refer to appropriate IFU for details.

To learn more about our Vital Signs Devices, please contact your local Hillrom representative or Customer Care Team, or europa@hillrom.com or middleeast@hillrom.com.



HillromTM

hillrom.com

The Welch Allyn® Vital signs devices are intended to be used by clinicians and medically qualified personnel for monitoring vital signs of all patients. Medical Device (93/42/EEC): Class IIa (Spot Vital Signs 4400 Device, Connex Spot Monitor, Connex Integrated Wall System) Class IIb (Connex Vital Signs Monitor), Notified Body: DQ. Notified Body: DQS Medizinprodukte GmbH (CE0297) Manufacturer: Welch Allyn, Inc., 4341 State Street Road, Skaneateles Falls, NY 13153 USA. EC REP: Welch Allyn Limited, Navan Business Park, Dublin Road, Navan, County Meath C15 AW22 Ireland. These medical devices are regulated health products which, pursuant to such regulation bears a CE mark. Carefully read the detailed instructions for safe and proper use included in the documents accompanying the medical devices. The personnel of healthcare establishments are responsible for the proper use and maintenance of these medical devices. Welch Allyn Inc. reserves the right to make changes without notice in design, specifications and models. The only warranty Welch Allyn Inc. makes is the express written warranty extended on the sale of its products.